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(EN) Instructions for use

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(SE) Bruksanvisning

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* 1 Mixingblock-Messeinheit = 1 cm (0,39 inch)

alphasil® PERFECT

Elastomeric precision impression material on polysiloxane basis, condensation-curing

Thank you for the confidence you have shown in our company. You have made the right decision with the alphasil® PERFECT impression system. Please be sure to read these instructions for use through carefully before starting to use the product!

alphasil® PERFECT products are only intended for use in dental applications by correspondingly qualified professionals. The medical device must be employed responsibly and with due caution.

Product description:

alphasil® PERFECT is a precise dental impression material on a C-silicone basis with an outstanding price-performance ratio which satisfies the highest demands. All alphasil® PERFECT products are made from carefully selected raw materials with state-of-the-art production technology in Germany.

Indication:

Crowns and bridges, partial crowns, inlays, onlays, complete and partial dentures.

Impression technique:

- **alphasil® PERFECT PUTTY SOFT:**
Basic material for the double mixing technique (1-step) and putty-wash technique (2-step)
- **alphasil® PERFECT HEAVY:**
Functional impressions, relining, basic material for double mixing technique
- **alphasil® PERFECT MEDIUM / LIGHT:** Wash material for putty-wash technique and double mixing technique

Dosage / processing: [see also Dosage table 2]

alphasil® PERFECT PUTTY SOFT

- **Dosage with alphasil® PERFECT ACTIVATOR PASTE:**
Remove the required quantity of putty and catalyst in the correct dosage ratio as shown in the illustration. Draw out the putty and catalyst to create a small strand, then fold it together again and knead it with your fingertips. Repeat this process a few times until you achieve a homogeneous colour, but always complying with the minimum mixing time of 30 seconds. Close containers and tubes again immediately after use!

alphasil® PERFECT HEAVY / MEDIUM / LIGHT

- **Dosage with alphasil® PERFECT ACTIVATOR PASTE:**
Place alphasil® PERFECT LIGHT or MEDIUM or HEAVY and activator paste side by side in the correct dosage ratio (2:1) as shown in the illustration on a mixing block or glass plate. Always ensure an even strand diameter. Pick up the activator with a spatula, distribute it evenly over the other component and mix thoroughly for 30 seconds until you achieve a homogeneous colour. Smooth the mixed material out with the spatula to remove any air.

Mixing ratio:

Adding too much or too little activator can affect the total setting time. As a general rule, adding more catalyst results in a faster setting time, whereas less catalyst slows it down. It is possible to increase or reduce the amount of activator by +/- 25%.

Impression procedure:

When using a suitable tray adhesive, please follow the manufacturer's instructions for use. After inserting the tray with the alphasil® PERFECT materials into the patient's mouth, press it into place briefly and then hold in place until fully cured. Always check the material intraorally to assess the degree of curing. When using the putty-wash technique (2-step), the first impression must be cleaned and dried in order to ensure cross-linking with the wash material.

Disinfection:

The impressions can be disinfected with a 2% glutaraldehyde solution. After removing the impression from the mouth, rinse it under running water for 15 seconds. Then disinfect the impression by immersing it in the disinfectant thoroughly. Observing the disinfectant manufacturer's specifications is essential. We recommend the (aldehyde-free) products Dentoprint MD liquid / Dentoprint MD pur from Müller-Omicron GmbH & Co. KG. Further information can be found at www.dgmk.de, and in J. Prosthet. Dent. 1999 May, 81(5), 621.

Instructions for use

Electroplating:

The impressions can be electroplated with copper or silver in standard baths.

Model fabrication:

The impression can be poured after just 30 mins and for up to 7 days. All class 3 and 4 gypsums which comply with the requirements in ISO 6873, e.g., Moldano, Kulzer (Type 3) and Fujirock EP, GC (Type 4), can be used. Undamaged impressions can be poured several times.

Working and safety instructions:

- The products described above must only be used in accordance with their instructions for use. The user has sole liability for any other use which contravenes these instructions for use.
- Close containers and tubes tightly immediately after use and do not interchange the lids and measuring spoons.
- Solutions used prior to impression taking, such as mouth rinses or retraction solutions, must be removed thoroughly.
- Do not leave any impression material residues in the oral cavity.
- Avoid contact with clothing, as cross-linked silicones are chemically resistant and leave stains which cannot be removed.
- Do not use alphasil® PERFECT products with polyether, polysulphides or addition-curing silicones.
- To prevent possible incompatibilities, do not combine alphasil® PERFECT products with C-silicones from other manufacturers.
- Storage: Store in a dry place at between 5°C and 27°C (41°F and 80.6°F).
- Shelf life: See date printed on package / label. Do not use after expiry date.
- Processing temperature: Room temperature between 18°C and 25°C (64.4°F to 77°F).

Contraindication:

The product must not be used in patients with known allergies to flavourings or other ingredients (see safety data sheet or ask the manufacturer). This should be clarified before use!

This product should not be used in pregnant women, breastfeeding mothers or children under six years of age. Take the impression with an A-silicone instead, e.g., betasil VARIO. This product may only be used by dental professionals.

Technical data and product features

Table 1:

Technical data	alphasil® PERFECT PUTTY SOFT	alphasil® PERFECT HEAVY	alphasil® PERFECT MEDIUM	alphasil® PERFECT LIGHT
Consistency / EN ISO 4823	Kneadable / Type 0	Heavy-bodied / Type 1	Medium-bodied / Type 2	Light-bodied / Type 3
Colour	White	Green	Blue	Dark blue
Mixing time	30 sec.	30 sec.	30 sec.	30 sec.
Working time incl. mixing time*	1 min. 45 sec.	1 min. 45 sec.	1 min. 45 sec.	1 min. 45 sec.
Time in mouth	≥ 4 min.	≥ 4 min.	≥ 4 min.	≥ 4 min.
Total setting time	5 min. 45 sec.	5 min. 45 sec.	5 min. 45 sec.	5 min. 45 sec.
Deformation under pressure (approx.)	2.0%	5.0%	6.0%	6.0%
Recovery from deformation (approx.)	98.8%	99.5%	99.7%	99.5%
Linear dimensional change (approx.)	0.20%	1.0%	1.0%	1.0%
Flavour	Mint	Mint	Mint	Mint

* At 23°C / 73.4°F. The times specified in the table above are based on a relative humidity of 50 +/- 10%.

As a general rule: higher temperatures accelerate setting and lower temperatures delay it.

Dosage table

Table 2:

	alphasil® perfect ACTIVATOR PASTE
alphasil® perfect PUTTY SOFT  1 ± 12 ml (22,4 g)	2,5*  ± 0,53 ml (0,5 g)
alphasil® perfect HEAVY  2* ± 1,22 ml (1,7 g)	1* 
alphasil® perfect MEDIUM  2* ± 1,16 ml (1,5 g)	1* 
alphasil® perfect LIGHT  2* ± 0,98 ml (1,18 g)	1*  ± 0,21 ml (0,2 g)

* 1 mixing block measuring unit = 1 cm [0.39 inches]

Warnings for alphasil® Activator Paste

- **Inhalation:** Take person out into the fresh air, administer oxygen if necessary.
- **Contact with skin:** Avoid prolonged contact with the skin. Following contact, wash the skin thoroughly with soap and water. We recommend wearing suitable protective clothing and gloves.

General warnings

- **Eye contact:** Avoid contact with eyes. In case of contact, rinse eyes immediately with plenty of clean water and consult an ophthalmologist without delay.
- **Swallowing:** If swallowed, rinse out the mouth and drink plenty of water. Always contact a doctor.
- **Allergic reactions:** It is not possible to rule out allergic reactions in hypersensitive people. If acute allergic reactions should occur during the treatment, treatment with the product should be stopped immediately and a doctor should be consulted. In case of known allergies, we generally recommend taking the impression with an A-silicone, e.g., betasil VARIO from Müller-Omicron GmbH & Co. KG.

Should a serious incident occur unexpectedly in connection with this medical device, it must be reported to the manufacturer and the competent authority in your Member State.

Please also observe the warnings in the safety data sheet!

Information at www.mueller-omicron.com

Errors and omissions excepted.

Contains hazardous substances:
diocetyl tin dilaurate

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